

tomo_nuevo

September 17, 2019

1 Test New Tomography Experiment Control

Test the OOP-based Tomography experiment control

```
[1]: import numpy as np
import matplotlib
import bluesky
import ophyd
import apstools
import databroker

from datetime import datetime
from seisidd.experiment import Tomography
```

```
[2]: # instantiate the experiment handle

testexp = Tomography()
```

/home/beams/S6BM/opt/miniconda3/envs/bluesky_py37_tornado6/lib/python3.7/site-packages/epics/pv.py:110: UserWarning: The `context` kwarg for epics.get_pv() is deprecated. New PVs will _not_ be created in the requested context.
'The `context` kwarg for epics.get_pv() is deprecated. New PVs '

```
[3]: testexp.mode = 'dryrun'
```

```
[4]: # expose internal RunEngine
RE = testexp.RE

# setup the Metadata as before
RE.md['beamline_id'] = 'APS 6-BM-A'
RE.md['versions'] = {}
RE.md['versions']['apstools'] = apstools.__version__
RE.md['versions']['bluesky'] = bluesky.__version__
RE.md['versions']['databroker'] = databroker.__version__
RE.md['versions']['matplotlib'] = matplotlib.__version__
RE.md['versions']['numpy'] = np.__version__
RE.md['versions']['ophyd'] = ophyd.__version__
RE.md['SESSION_STARTED'] = datetime.isoformat(datetime.now(), " ")
```

```
[5]: scan_cfg = 'seisidd/config/tomo_scan_template.yml'
```

```
[6]: # summarize plan  
testexp.dryrun(scan_cfg)
```

===== Open Run =====

```
A_shutter -> open  
det_tiff1_file_path -> /dev/shm/tmp/  
det_tiff1_file_name -> ttt  
det_tiff1_file_write_mode -> 2  
det_tiff1_num_capture -> 26  
det_tiff1_file_template -> %s%s_%06d.hdf  
det_hdf1_file_path -> /dev/shm/tmp/  
det_hdf1_file_name -> ttt  
det_hdf1_file_write_mode -> 2  
det_hdf1_num_capture -> 26  
det_hdf1_file_template -> %s%s_%06d.hdf  
det_tiff1_enable -> 0  
det_hdf1_enable -> 1  
det_hdf1_capture -> 1  
det_cam_frame_type -> 0  
tomostage_ksamX -> -1.0  
tomostage_ksamZ -> 0.0  
det_hdf1_nd_array_port -> PROC1  
det_tiff1_nd_array_port -> PROC1  
det_proc1_enable -> 1  
det_proc1_reset_filter -> 1  
det_proc1_num_filter -> 1  
det_cam_trigger_mode -> Internal  
det_cam_image_mode -> Multiple  
det_cam_num_images -> 5  
  Read ['det']  
tomostage_ksamX -> 0.0  
tomostage_ksamZ -> 0.0  
det_cam_frame_type -> 1  
det_hdf1_nd_array_port -> PROC1  
det_tiff1_nd_array_port -> PROC1  
det_proc1_enable -> 1  
det_proc1_reset_filter -> 1  
det_proc1_num_filter -> 1  
tomostage_preci -> 0.0  
  Read ['det']  
tomostage_preci -> 0.5  
  Read ['det']  
tomostage_preci -> 1.0  
  Read ['det']  
tomostage_preci -> 1.5  
  Read ['det']
```

```

tomostage_preci -> 2.0
  Read ['det']
tomostage_preci -> 2.5
  Read ['det']
tomostage_preci -> 3.0
  Read ['det']
tomostage_preci -> 3.5
  Read ['det']
tomostage_preci -> 4.0
  Read ['det']
tomostage_preci -> 4.5
  Read ['det']
tomostage_preci -> 5.0
  Read ['det']
det_cam_frame_type -> 2
tomostage_ksamX -> -0.9961946980917455
tomostage_ksamZ -> 0.08715574274765817
det_hdf1_nd_array_port -> PROC1
det_tiff1_nd_array_port -> PROC1
det_procl_enable -> 1
det_procl_reset_filter -> 1
det_procl_num_filter -> 1
det_cam_trigger_mode -> Internal
det_cam_image_mode -> Multiple
det_cam_num_images -> 5
  Read ['det']
tomostage_ksamX -> 0.0
tomostage_ksamZ -> 0.0
det_cam_frame_type -> 3
A_shutter -> close
det_hdf1_nd_array_port -> PROC1
det_tiff1_nd_array_port -> PROC1
det_procl_enable -> 1
det_procl_reset_filter -> 1
det_procl_num_filter -> 1
det_cam_trigger_mode -> Internal
det_cam_image_mode -> Multiple
det_cam_num_images -> 5
  Read ['det']
===== Close Run =====

```

[8]: `testexp.run(scan_cfg)`

```

Transient Scan ID: 2      Time: 2019-08-22 12:47:10
Persistent Unique Scan ID: 'e4c9c153-1090-47f1-bf23-fc246ac11792'
New stream: 'primary'
+-----+-----+

```

seq_num	time
1	12:47:16.2
2	12:47:22.3
3	12:47:25.4
4	12:47:28.4
5	12:47:31.5
6	12:47:34.6
7	12:47:37.6
8	12:47:40.7
9	12:47:43.7
10	12:47:46.8
11	12:47:49.9
12	12:47:52.9
13	12:47:57.5
14	12:48:02.1

generator tomo_scan ['e4c9c153'] (scan num: 2)

[8]: ('e4c9c153-1090-47f1-bf23-fc246ac11792',)

[]: